

# Barron Bala Asymmetric Cornea Punch



## Asymmetric punch shape:

- Novel sculpted blade design allows blade to fit curvature of the punch and allows the pedicle/tail to be cut last
- Notch/bump: Assists in orientation of the graft in the eye

## Features an “open top” design:

- Improves tissue visualization prior to cutting

## Optimal intended use:

- Surgeon self-preparation of donor corneal tissue
- Eye Bank preparation of donor corneal tissue

## Item #: K20-2206

- Punch diameter: 7.6mm

Barron Bala Asymmetric  
Cornea Vacuum Punch

