


Relieving the Pressure of Subretinal Injections!

PolyVent™ Cannulas

- Multiple egress vents allow pressure release, minimizing pressure on the bleb during subretinal injections.
- Green 5 mm long, 38 gauge tip for better visualization.
- 100% flow tested during manufacturing for reliable performance.
- Available in 23 gauge and 25 gauge.



3305 PolyVent™ Cannula 25g/38g
3306 PolyVent™ Cannula 23g/38g



- Available with extendable 38 gauge tip for easier insertion through valved cannulas.
- Small bore infusion tubing significantly reduces dead space to 0.1 ml, decreasing waste of injectable by 50% versus comparable subretinal devices.



3307 Extendable PolyVent™ Cannula 25g/38g
3308 Extendable PolyVent™ Cannula 23g/38g

 **MedOne**
670 Tallevast Road
Sarasota, FL 34243 USA
Tel: 941.359.3129
www.MedOne.com

©2021 MedOne Surgical, Inc (657)