

# assistPro

assistPro F Item No.: 2400 0005  
 assistPro H Item No.: 2400 0009

The perfect assistant's chair for your OR.

## Chassis

- High-grade 5-arm aluminum pedestal
- Excellent handling characteristics thanks to smooth-running double castors
- Optionally available with braking double castors
- Supports loads up to 130 kilograms



## Seating system

### Comfort Seat

- Ergonomically shaped seat and backrest
- Height and angle adjustable lumbar support
- Supports a correct, back-friendly sitting posture



- Durable, easy-care seat cover made of medical-grade artificial leather to meet the highest OR requirements
- Optionally available with electrically conductive upholstery

**assistPro F** foot trigger  
sterile height adjustment



**assistPro H** manual trigger  
classic height adjustment



## Technical specifications

### Lifting column

- Height adjustment *assistPro F* 530 – 730 mm
- assistPro H* 520 – 720 mm

### Chassis

- Five legged chassis  $\varnothing$  585 mm
- Castors  $\varnothing$  65 mm

### Comfort seat

- D x W x H (incl. back part) 385 x 425 x 540 mm
- Seat angle adjustment +5°/-5°
- Back angle adjustment +5°/-25°

### Weight

- assistPro F* 12 kg
- assistPro H* 11 kg

### Max. load

130 kg

Subject to technical modifications.

## Flexibility

### Variants

- Sterile, comfortable height adjustment via foottriggered gas pressure spring (*assistPro F*)
- Flexible height adjustment from 53 to 73 cm (*assistPro F*)
- Classic seat height adjustment via manual trigger underneath the seat (*assistPro H*)
- Flexible height adjustment from 52 to 72 cm (*assistPro H*)
- Reduces neck and back tension caused by prolonged sitting
- Optionally available with arm supports (*assistPro F*)

### Standard colour

schwarz  
 Item No.: 3037 0235



**Comfort**

"Periline"

Page 111

# Components



*assistTrend*  
Item No.: 2500 0002



*assistPro F*  
Item No.: 2400 0005



*assistPro H*  
Item No.: 2400 0009



Armrest pair  
incl. traverse  
*assistTrend*

Item No.: 3300 0000



Armrest pair  
incl. traverse  
*assistPro F*

Item No.: 3351 0038



*easyPro F*  
Item No.: 2400 0010



*easyPro H*  
Item No.: 2400 0006



*easyPro R*  
Item No.: 2400 0013



*assistLine*  
Item No.: 2400 0015



Double rimmed castor,  
65 mm, with braking function

Item No.: 3351 0040

# "Periline"



## Product Specification

- Without Phthalates
- Biocompatibility:  
Cytotoxicity (EN ISO 10 993-5) and  
test for irritation and delayed-type  
hypersensitivity (EN ISO 10 993-10)
- Urine resistance
- Blood resistance
- Resistance to perspiration
- Resistance to spitting
- Light fastness
- Color fastness to rubbing
- Flame retardant

<i>surgiForce</i>	Comfort seat	2600 0004
<i>surgiForce light</i>	Comfort seat	2600 0006
<i>surgiTrend</i>	Comfort seat	2200 0012
<i>surgiLine</i>	Comfort seat	2200 0013
<i>assistPro F</i>		2400 0005
<i>assistPro H</i>		2400 0009
<i>easyPro F</i>		2400 0010
<i>easyPro H</i>		2400 0006
<i>easyPro R</i>		2400 0013

+ Item No.: 3037 0...



kastanie 015817  
Item No.: .....240



light brown 016037  
Item No.: .....239



creme 015812  
Item No.: .....225



weiß 010881  
Item No.: .....224



stein 015813  
Item No.: .....236



saphir 015814  
Item No.: .....231



delft 010766  
Item No.: .....241



antistatisch schwarz  
Item No.: .....223



schwarz 011499  
Item No.: .....235



anthracit 016046  
Item No.: .....234



azur 015815  
Item No.: .....230



aubergine 016041  
Item No.: .....229



lavendel 016042  
Item No.: .....228



rot 010878  
Item No.: .....238



orange 016040  
Item No.: .....237



kiwi 015816  
Item No.: .....233



lime green 016043  
Item No.: .....232



lemon 016039  
Item No.: .....226



safran 010879  
Item No.: .....227

The colours displayed on your device may differ slightly from the original colours. If you wish to receive original colour charts, please contact us at [info@ufsk-osys.com](mailto:info@ufsk-osys.com).

# Upgrade your OR to perfection.



*eyeForce*



*500 XLE*



*500 XLE comfort*



*500 bariatric*



*500 XLS*



*500 ECO*



*400 XL*



*motionLine*



*motionLine-E*



*surgiForce*



*surgiForce light*



*surgiTrend*



*surgiLine*



*stayClean wipes*



*stayClean*



*staySoft*



*stayClean roll*



*assistTrend*



*assistPro F*



*assistPro H*



*easyPro F*



*easyPro H*



*easyPro R*



*assistLine*



*supportLine*



*supportLine plus*

## CONTACT US



**UFSK**  
INTERNATIONAL

**OSYS**

UFSK-International OSYS GmbH  
Kirchhoffstraße 1  
93055 Regensburg  
Germany

Phone +49 (0)941 78862 15  
Fax +49 (0)941 78862 35  
Mail [info@ufsk-osys.com](mailto:info@ufsk-osys.com)  
Web [ufsk-osys.com](http://ufsk-osys.com)

**MADE IN GERMANY**  
Design | Construction | Manufacturing