

# assistTrend

Item No.: 2500 0002

The ergonomic assistant chair.  
The innovation for a healthy back.

## Chassis

- High-grade 5-arm aluminium pedestal
- Excellent handling characteristics thanks to smooth-running double castors
- Optionally available with braking double castors
- Supports loads up to 130 kilograms

## Seating system

### Dynamic seat

- Innovative seating system Dynamic for a healthy lower back
- Height and angle adjustable lumbar support
- Supports a correct, back-friendly sitting posture
- Durable, easy-care seat cover made of medical-grade artificial leather to meet the highest OR requirements



The seating system Dynamic has a lumbar support which is activated automatically by the surgeon's bodyweight. The intensity of ergonomic support can be adjusted underneath the seat.



## Technical specifications

### Lifting column

- Height adjustment 580 – 780 mm
- For body height of 1.75 – 1.95 m

### Chassis

- Five legged chassis ø 585 mm
- Castors ø 65 mm

### Dynamic seat

- D x W x H (incl. back part) 390 x 440 x 410 mm
- Seat/back angle adjustment +12°

### Weight

7.5 kg

### Max. load

130 kg

Subject to technical modifications.

## Flexibility

- Sterile, comfortable height adjustment via foottriggered gas pressure spring
- Flexible height adjustment from 58 to 78 cm
- Reduces neck and back tension caused by prolonged sitting
- Optionally available with arm supports
- Helps maintain sterile conditions

### Standard colour

schwarz

Item No.: 3037 0021



**Dynamic**

"Stamskin Top"

Page 110

# Components



*assistTrend*  
Item No.: 2500 0002



*assistPro F*  
Item No.: 2400 0005



*assistPro H*  
Item No.: 2400 0009



Armrest pair  
incl. traverse  
*assistTrend*

Item No.: 3300 0000



Armrest pair  
incl. traverse  
*assistPro F*

Item No.: 3351 0038



*easyPro F*  
Item No.: 2400 0010



*easyPro H*  
Item No.: 2400 0006



*easyPro R*  
Item No.: 2400 0013



*assistLine*  
Item No.: 2400 0015



Double rimmed castor,  
65 mm, with braking function

Item No.: 3351 0040

# "Stamskin Top"



User friendly: High flexibility of the material ensures a perfect finish.

Safe for the user: Phthalate-free, easy to clean and fulfills fire safety standards.

## Product Specification

- Extremely durable
- Phthalate free
- Easy to clean
- Flame protection
- Mildew resistant

*surgi*Trend  
*assist*Trend

Dynamic seat 2200 0010  
2500 0002

+ Item No.: 3037 0...



F4340-10120  
Item No.: .....023



F4340-07422  
Item No.: .....016



F4340-20130  
Item No.: .....084




F4340-20156  
Item No.: .....207



F4340-50622  
Item No.: .....215




F4340-50623  
Item No.: .....216




F4340-20129  
Item No.: .....146




F4340-00002  
Item No.: .....021



F4340-20162  
Item No.: .....205



F4340-20161  
Item No.: .....206



F4340-20150  
Item No.: .....204



F4340-5015  
Item No.: .....198



F4340-07496  
Item No.: .....046



F4340-07432  
Item No.: .....002



F4340-20154  
Item No.: .....208



F4340-07445  
Item No.: .....030



F4340-50621  
Item No.: .....199



F4340-07480  
Item No.: .....039



F4340-07436  
Item No.: .....020



F4340-10295  
Item No.: .....044



F4340-07459  
Item No.: .....028



F4340-20157  
Item No.: .....211



F4340-5347  
Item No.: .....209



F4340-50454  
Item No.: .....210



F4340-20149  
Item No.: .....212



F4340-50619  
Item No.: .....214



F4340-20127  
Item No.: .....079



F4340-07477  
Item No.: .....036



F4340-07476  
Item No.: .....035



F4340-50618  
Item No.: .....220



F4340-07485  
Item No.: .....034



F4340-50620  
Item No.: .....213



F4340-07478  
Item No.: .....037



F4340-07479  
Item No.: .....038

# Upgrade your OR to perfection.



*eyeForce*



*500 XLE*



*500 XLE comfort*



*500 bariatric*



*500 XLS*



*500 ECO*



*400 XL*



*motionLine*



*motionLine-E*



*surgiForce*



*surgiForce light*



*surgiTrend*



*surgiLine*



*stayClean wipes*



*stayClean*



*staySoft*



*stayClean roll*



*assistTrend*



*assistPro F*



*assistPro H*



*easyPro F*



*easyPro H*



*easyPro R*



*assistLine*



*supportLine*



*supportLine plus*

## CONTACT US



**UFSK**  
INTERNATIONAL

**OSYS**

UFSK-International OSYS GmbH  
Kirchhoffstraße 1  
93055 Regensburg  
Germany

Phone +49 (0)941 78862 15  
Fax +49 (0)941 78862 35  
Mail [info@ufsk-osys.com](mailto:info@ufsk-osys.com)  
Web [ufsk-osys.com](http://ufsk-osys.com)

**MADE IN GERMANY**  
Design | Construction | Manufacturing