

assistPro

assistPro F Item No. : 2400 0005
 assistPro H Item No. : 2400 0009

The perfect assistant's chair for your OR.

Chassis

- High-grade 5-arm aluminum pedestal
- Excellent handling characteristics thanks to smooth-running double castors
- Optionally available with braking double castors
- Supports loads up to 130 kilograms



Seating system

Comfort Seat

- Ergonomically shaped seat and backrest
- Height and angle adjustable lumbar support
- Supports a correct, back-friendly sitting posture



- Durable, easy-care seat cover made of medical-grade artificial leather to meet the highest OR requirements
- Optionally available with electrically conductive upholstery



assistPro F foot trigger
sterile height adjustment



assistPro H manual trigger
classic height adjustment



Technical specifications

Lifting column

Height adjustment	
- assistPro F	530 - 730 mm
- assistPro H	520 - 720 mm

Chassis

- Five legged chassis	ø 585 mm
- Castors	ø 65 mm

Comfort seat

- D x W x H incl. backrest	385 x 425 x 540 mm
- Seat angle adjustment	+5°/-5°
- Back angle adjustment	+5°/-25°

Weight

assistPro F	12 kg
assistPro H	11 kg

Max.load

130 kg

Subject to technical modifications.

Flexibility

Variants

- Sterile, comfortable height adjustment via foottriggered gas pressure spring (*assistPro F*)
- Classic seat height adjustment via manual trigger underneath the seat (*assistPro H*)
- Flexible height adjustment from 52 to 71.5 cm
- Reduces neck and back tension caused by prolonged sitting
- Optionally available with arm supports (*assistPro F*)



Comfort

"Periline"
Page 111

Accessories



assistTrend
Item No. : 2500 0002



assistPro F
Item No. : 2400 0005



assistPro H
Item No. : 2400 0009



assistTrend
Arm rest pair
incl. traverse

Item No.: 3300 0000



assistPro F
Arm rest pair
incl. traverse

Item No.: 3351 0038



easyPro F
Item No. : 2400 0010



easyPro H
Item No. : 2400 0006



easyPro R
Item No. : 2400 0013



assistLine
Item No. : 2400 0015



Double rimmed caster
65 mm
with braking function

Item No.: 3351 0040

"Periline"



surgiForce
Comfort seat 2600 0004

surgiForce light
Comfort seat 2600 0006

surgiTrend
Comfort seat 2200 0012

surgiLine
Comfort seat 2200 0013

assistPro F 2400 0005

assistPro H 2400 0009

easyPro F 2400 0010

easyPro H 2400 0006

easyPro R 2400 0013



kastanie 015817
Item No.:240



light brown 016037
Item No.:239



creme 015812
Item No.:225



weiß 010881
Item No.:224



stein 015813
Item No.:236



saphir 015814
Item No.:231



delft 010766
Item No.:241



antistatisch schwarz
Item No.:223



schwarz 011499
Item No.:235



anthracit 016046
Item No.:234



azur 015815
Item No.:230



aubergine 016041
Item No.:229



lavendel 016042
Item No.:228



rot 010878
Item No.:238



orange 016040
Item No.:237



kiwi 015816
Item No.:233



lime green 016043
Item No.:232



lemon 016039
Item No.:226



safran 010879
Item No.:227

+ Item No.: 30370....

The printed colours may differ slightly from the original colours. If you wish to receive original colour charts, please contact us at info@ufsk-osys.com.

Upgrade your OR to perfection.



eyeForce



500 XLE



500 XLE comfort



500 bariatric



500 XLS



500 ECO



400 XL



motionLine



motionLine-E



surgiForce



surgiForce light



surgiTrend



surgiLine



stayClean wipes



stayClean



staySoft



stayClean roll



assistTrend



assistPro F



assistPro H



easyPro F



easyPro H



easyPro R



assistLine



supportLine



supportLine plus

CONTACT US



UFSK-International OSYS GmbH
Kirchhoffstraße 1
93055 Regensburg
Germany

Phone +49(0)941 78862 15
Fax +49(0)941 78862 35
Mail info@ufsk-osys.com
Web ufsk-osys.com

MADE IN GERMANY
Design | Construction | Manufacturing