



Visitec® Single Use Manipulators (Steel Handle)

Single-Use Sterile Material

Distribuido por:



Product class

Class IIa, sterile, single-use
According to European Medical Device
Directive 93/42/EEC

Manufacturer

Beaver-Visitec International Ltd
Bidford, UK

Certification Information

•CE Certificate

CE 0120
SGS Certificate US96/7138

•ISO Certificates

ISO 9001:2008 - Certificate US96/7137
ISO 13485:2003 - Certificate US96/7140
ISO 14001:2004 - Certificate 423

Sterile product

- Gamma irradiation (cobalt 60)
- Sterigenics, Belgium

Shelf life

- 5 years from manufacturing date

Indications

Single-use instruments for cataract surgery,
designed to manipulate the nucleus or the
intraocular lens in the capsular bag.

Product composition

- Tip:** Stainless Steel
- Handle:** Stainless Steel
- Absence of:** Latex / PVC / Phthalates/
Animal origin components

Packaging material

- Blister:** Coated Tyvek, PETG
- Shipper:** Paper carton








General features

- Made of stainless steel: precision with the convenience of single use.
- Visitec® Single Use Instrument help reduce the risk of infection transmission and help you comply with today's safety equipment's and recommendation.
- Easy to use 100% traceably, eliminate the need for costly of sterilization

Warning & Precautions: See box label and insert.

Ordering Information:

Visitec® Single Use Manipulators (Steel Handle)		
Item #	Description	Picture
BVI-581412	Burrato's Hook, angled, 0.25 mm diameter, blunt tip. Overall length 120 mm. Use to manipulate the nucleus.	
BVI-581423	Micro Manipulator, angled, 0.35 mm diameter, blunt tip. Overall length 120 mm. Used to manipulate the nucleus.	
BVI-581464	Lens Manipulator [Lester], angled, 0.44 mm diameter hourglass shape tip. Overall length 115 mm. Use for IOL manipulation.	
BVI-581465	Nucleus Manipulator [Drysdale], flattened paddle tip, 1.10 mm tip. Overall length 115 mm. Used to manipulate the nucleus.	
BVI-581466	Nucleus Rotator [Pisacano], Angled, 0.70 mm diameter flat, blunt tip. Overall length 115 mm. Used to manipulate the nucleus.	

Distribuido por:

MEDICALMIX

PASIÓN POR LA OFTALMOLOGÍA

Amposta, 20 | 08174 Sant Cugat del Vallès
T. 935 893 637 | F. 935 892 829
mmix@medicalmix.com | www.medicalmix.com


Keeping Your Vision in Sight

24 rue Lamartine
F-38320 Eybens
Tel.: 0033 4 85 88 00 33
Fax: 0033 4 56 80 92 27
www.beaver-visitec.com